Operation Manual

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BeWell®

BeWell Device (Hand-held)

Introduction

BeWell Laser Device, is an innovative, high-tech health care product developed and introduced from Germany. The device releases 650 nanometer low level laser, which is easy to carry and apply. It offers therapy that is non-invasive, painless, non-medicinal and with no side-effects. Regular use of this device will improve red blood cell function and microcirculation. It acts by instantly increasing erythrocyte deformability, inhibiting the aggregation of erythrocytes and platelets, reducing blood viscosity and increasing the oxygen-transporting capacity of the red blood cells. Moreover, low level laser has been shown to promote the process of anti-oxidation and improve immune function.

Using BeWell Device daily can help effectively control blood pressure and alleviate the side-effects of lipid-lowering drugs. It also enhances the activity of insulin receptors and improves insulin function in diabetic patients. Not only can it restore blood's beta-endorphin levels, it can also completely prevent symptoms of premenstrual syndrome. It is especially effective in preventing and improving conditions like coronary heart disease, stroke, hypertension, hyperlipidemia, diabetes, premenstrual syndrome (PMS) and other common ailments. In addition, this device is effective in relieving and preventing various sub-health symptoms which include headache, dizziness, tinnitus, insomnia, nasal allergies, dry throat, skin allergies, back pain, anxiety, depression, and other mood disorders.

Thank you for purchasing our product and welcome to BeWell. Please read this operation manual carefully before you start to use the device. This manual shows the detail operation and user guide.

 
BeWell[®]

Table of Contents

1. Important Safety Instructions
2. Product Layout
3. Package Content
4. User Guide
5. Care and Maintenance
6. Product Specifications
7. Service and Warranty

Important Safety Instructions

1

- 1.1 Do not attempt to open the device unless otherwise specified in this operation manual.
- 1.2 Do not disassemble or drop the device.
- 1.3 Keep the device away from high temperature and vibration, avoid exposing it under direct sunlight.
- 1.4 Do not expose the device to rain or any type of moisture.
- 1.5 Avoid exposure to laser beam other than indicated.
- 1.6 Keep the device and all its parts including accessories out of the reach of small children.
- 1.7 If you have any operation problem, please contact our Customer Service Team.
- 1.8 Battery needs to be charged for 8 hours before first time of operation.

Type BF applied part

Attention, consult accompanying documents

Disposal — Do not dispose the electrical appliances as unsorted municipal waste. Use separate collection facilities. Contact your local government for information regarding the collection systems available. If electrical appliances are disposed in landfills or dumps, hazardous substances can leak into the groundwater and get into the food chain, it can damage your health and well-being.

People with **blood coagulation dysfunction** (e.g. vascular purpura, thrombocytopenic purpura, haemophiliac), **pregnancy, menstruation** and **electrical pacemaker of the heart** should use carefully, and consult your doctor before using the device.

Product Layout

2

Package Content

3

3.1 Cradle

3.2 Adaptor Set

3.3 Operation manual

3.4 Travel bag for the device

User Guide

4

A. Get Started

- 4.1 Hold the device in your left hand, make sure the emitter has close contact at the web space area of your hand.

- 4.2 Turn on the device by pressing the  button for 2 seconds, the therapy time indicator and battery level indicator will turn on in blue.

- 4.3 Confirm and start by pressing OK button, therapy time indicator will flash slowly.
- 4.4 The treatment starts when laser is turned on.

- 4.5 After 30-minutes treatment, a '3-beep' alarm will sound and the device will switch off automatically.
- 4.6 If the device is not held in the right position, a 'beep' alarm will sound and the device will switch off automatically.

B. Indicating the treatment time

- 4.7 Each treatment cycle lasts for 30 minutes.
- 4.8 Display the treatment remaining time by pressing the OK button while the device is operating.
- 4.9 Each LED sector of the Therapy Time Indicator represents the remaining treatment time.

30 minutes left

20 minutes left

10 minutes left

30 minutes treatment is done

C. Turning off the device

- 4.10 Turn off the device by pressing the  button for 3 seconds, the device will switch off automatically.
- 4.11 The device will switch off automatically when it is not held in the right position in your hand.

D. Indicating the battery level

- 4.12 Each LED sector of the battery level indicator shows the remaining battery level.

Fully Charged

4/5 left

3/5 left

2/5 left

1/5 left

Please recharge the battery as soon as possible when the LED Sector turns orange.

E. Recharging battery

4.13 Connect the cradle to the adaptor and put the device into the cradle.

4.14 To indicate the battery is charging, the battery level indicator flashes with blue light and the LED light at the cradle lights up in orange.

4.15 After the battery is fully charged, the battery level indicator will turn off, while the orange LED light at the cradle will turn blue.

Note: The battery can also be charged by using provided adaptor without cradle

Note: To keep the battery fully charged, you may keep the device at the cradle when not in use.

Warning: Do not operate the device when it is charging.

F. Silent Mode

- 4.16** Silent Mode must be set while the device is operating.
- 4.17** Mute : Press the OK button until heard a 'Beep' sound once; NO alarm will be sound after 30-minutes treatment and if the device is not held in the right position.
- 4.18** Vocal : Press the OK button until a 'Beep' sound is heard twice; alarm will be sound after 30-minutes treatment and if the device is not held in the right position.

Care and Maintenance

5

- 5.1** Keep the device away from precipitation. Humidity and liquids contain minerals that will corrode electronic circuits.
- 5.2** The device is a product of superior design and craftsmanship and should be treated with care. Proper operation allows you to enjoy this product longer. For any question about operation, please contact our customer service team. We are glad to offer our assistance and professional services.

Product Specifications

6

Name of Product	: BeWell Device (Hand-held)
Laser wavelength	: 650nm
Maximum laser output	: 5 mW x 3
Preset laser output	: 5 mW
Battery	: Lithium
Power Consumption	: 3.7VDC 720mA
Power Supply	: Input 100V – 240VAC 50/60Hz Output 5.5VDC 1A
Safety Classification	: Electrical appliance safety type II
Safe mark	:  Laser radiation, avoid exposure to beam, class 3R laser product

- 7.1 If the device does not work or does not work properly, please check the list below

Problem

- The laser does not turn on
- The device turns off automatically

Possible Cause(s)

- The "OK" button has not been pressed to activate the device
- The device has not been held in hand properly

Solution(s)

- Press "OK" button to activate the laser
- Check if the device is held in a right position

- 7.2 If the problem is not mentioned in the list or still not working, please contact our Customer Service Center.
- 7.3 This product consists of intricate electrical parts, do not attempt to open the device. If the device is not working properly, take it to our Customer Service Center.
- 7.4 BeWell Global Limited warrants each new hand-held device against defects in materials or workmanship for a period of two years from the date of purchase, and agrees to repair or replace any defective product.

IMPORTANT:

This warranty does not cover damage resulting from accident, misuse or abuse, lack of reasonable care and unauthorized disassembly.

Copyright

This manual is copyright by BeWell, and its licensee. Under the copyright laws, this manual may not be reproduced in any form, in whole or part without the prior consent of BeWell and its licensee.

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